



## Z-Band UHF RFID Wristbands

### Improve Accuracy and Operational Efficiency

Zebra's Z-Band® UHF RFID durable direct thermal printable wristband solutions simplify the identification of patients in clinical environments or guests at event venues. These disposable RFID wristbands are specifically designed for optimal use with Zebra's ZD500R desktop printer and Zebra RFID readers. Zebra's UHF RFID wristbands offer a wristbanding solution to fit your RFID application, helping you better manage events or ensure positive patient identification.

### Z-Band Direct RFID SR Wristband

Zebra's Z-Band Direct RFID SR wristband is best for use cases where the wristband will be in close proximity to the RFID reader.

### Common Applications

#### Healthcare\*

- At the patient bedside through a blanket
- During a patient's surgery through the drape
- Allows caregivers to obtain immediate access to critical patient records and update them at the bedside
- Quickly associate supplies and equipment used in treatments

#### Hospitality

- Guest Identification
- Access Management
- Cashless POS



### What Can Passive UHF RFID Do?

Passive UHF RFID technology can utilize short and long range RFID readers to offer the tracking of any tagged asset with last known read location throughout a facility.

| Part Number | Color | Size     | Inlay  | Read Range | Latex-free |
|-------------|-------|----------|--------|------------|------------|
| 3014578     | White | 1" x 11" | BT0600 | Up to 3m   | Yes        |

**FACT SHEET**

**Z-BAND UHF RFID WRISTBAND**

**Z-Band UltraSoft RFID LR Wristband**

Zebra’s Z-Band UltraSoft LR Wristband is uniquely designed with a specialty inlay designed for use against the human body. The design allows the wristband to flag away from the body, enabling longer read ranges. Ideal for use when overhead readers or portals are being used to provide best last known location tracking abilities.



**Common Applications**

**Healthcare\***

- Streamlining patient tracking throughout a medical facility

**Hospitality**

- Guest tracking throughout a venue

| Part Number  | Color | Size             | Inlay   | Read Range   | Latex-free |
|--------------|-------|------------------|---------|--------------|------------|
| ZIPRD3015155 | White | 1.1875” x 14.25” | ZBR2002 | Up to 17m ** | Yes        |

**Custom Solutions**

Zebra can create a custom wristband featuring your preferred inlay and chip combination, as well as pre-printed logos and colors.

\*Wristband should be removed from the patient during an MRI.

\*\*Read range is highly contingent on where inlay is placed in relation to the body.

**For more information, visit [www.zebra.com/wristbands](http://www.zebra.com/wristbands)**



**NA and Corporate Headquarters**  
+1 800 423 0442  
inquiry4@zebra.com

**Asia-Pacific Headquarters**  
+65 6858 0722  
contact.apac@zebra.com

**EMEA Headquarters**  
zebra.com/locations  
contact.emea@zebra.com

**Latin America Headquarters**  
+1 786 245 3934  
la.contactme@zebra.com